

## RIH – STEREOTACTIC GAMMA KNIFE GE LIGHTSPEED 16 / OPTIMA CT580 PROTOCOL

**Indications:** A gamma knife is a neurosurgical device used to treat brain tumors with radiation therapy.

<b>Position/Landmark</b>	Supine head first, in stereotactic headholder. Neuro team will position.				
<b>Topogram Direction</b>	Craniocaudal				
<b>Respiratory Phase</b>	Any				
<b>Scan Type</b>	Helical				
<b>KV / mA / Rotation time (sec)</b>	120kv / smart mA (50-250) / 0.8 sec				
<b>Pitch / Speed (mm/rotation)</b>	.562:1 , 5.62mm				
<b>Noise Index</b>	5.00				
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 16 = 10mm				
<b>Average Tube Output</b>	ctdi – 51.1 mGy dlp – 872 mGy.cm				
<b>Helical Set</b>					recon
Slice Thickness/ Spacing	recon	body part	thickness/ spacing	algorithm	recon destination .
Algorithm	1	thin gamma knife	0.6mm x 0.6mm	bone	gamma knife
Recon Destination	2	brain	5mm x 5mm	standard	pacs
	3	skull	5mm x 5mm	bone	pacs
<b>Scan Start / End Locations</b>	Prescribed by gamma knife.				
<b>DFOV</b>	25cm decrease appropriately				
<b>IV Contrast Volume / Type / Rate</b>	The gamma knife team will prescribe iv contrast only if necessary. Typical iv contrast dosage for this protocol is 200cc omni 350, hand injected.				
<b>Scan Delay</b>	Prescribed by gamma knife.				
<b>2D/3D Technique Used</b>					
<b>Comments:</b> Recon 1 is a bone algorithm, thin data set sent to gamma knife. Recon 2 is a standard algorithm brain, Recon 3 is for the skull.					
<b>Images required in PACS</b>	Scouts, 5mmx 5mm standard brain, 5mm x 5mm skull, Dose Report				