

**RIH – IV CONTRAST NECK  
SIEMENS DEFINITION AS+ PROTOCOL**

**Indications - mass, lymphoma, adenopathy, mets.**

<b>Position/Landmark</b>	Head first or feet first-Supine 1cm superior to skull vertex			
<b>Topogram Direction</b>	Craniocaudal / Craniocaudal			
<b>Respiratory Phase</b>	Inspiration			
<b>Scan Type</b>	Helical			
<b>Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization</b>	Care kV 120 / Care Dose4D 115 / 1 sec .8:1 , 32.00mm 3 / 7			
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 64 = 40mm (128 x .6mm)			
<b>Average Tube Output</b>	ctdi – 9.5 mGy dlp – 300 mGy.cm			
<b>Helical Set</b>	body	thickness/		recon
Slice Thickness/ Spacing	recon	part	spacing	algorithm
Algorithm	1	<b>axial iv neck</b>	3mm x 3mm	I40f medium
Recon Destination	2	<b>coronal iv neck</b>	3mm x 3mm	I40f medium
	3	thin neck	.75mm x .7mm	I40f medium
<b>Scan Start / End Locations</b>	External auditory meatus Aortic arch			
<b>DFOV</b>	20 cm decrease appropriately			
<b>IV Contrast Volume / Type / Rate</b>	70mL Iohexol (Omnipaque 300) , 2mL/sec if needed			
<b>Scan Delay</b>	35 seconds			
<b>2D/3D Technique Used</b>	Workstream 4D mpr of 3mm x 3mm <b>coronal neck</b> series, auto-transferred to PACS.			
<b>Comments:</b> Recon 3 is a thin helical volume of the neck that is archived to the TeraRecon server.				
<b>Images required in PACS</b>	Topograms, 3mm x 3mm neck, 3mm x 3mm coronal neck, Patient Protocol			