

**RIH – LUMBOSACRAL SPINE
SIEMENS DEFINITION AS+ PROTOCOL**

Indication: fracture, trauma, mets, disc rupture, disc herniation, stenosis, post myelogram.

Position/Landmark	Head first or feet first-Supine Mid Skull																														
Topogram Direction	Craniocaudal / Craniocaudal																														
Respiratory Phase	Suspension																														
Scan Type	Helical																														
Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization	Care kV 120 / Care Dose4D 230 / 1.0 sec .8:1 , 32.00mm 2 / 3																														
Detector width x Rows = Beam Collimation	0.625mm x 64 = 40mm (128 x .6mm)																														
Average Tube Output	ctdi – 10.0mGy dlp – 300mGy.cm																														
Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>axial l spine tissue</td> <td>3mm x 3mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>2</td> <td>axial l spine bone</td> <td>3mm x 3mm</td> <td>I70h very sharp</td> <td>pac</td> </tr> <tr> <td>3</td> <td>coronal l spine</td> <td>3mm x 3mm</td> <td>I70h very sharp</td> <td>pac</td> </tr> <tr> <td>4</td> <td>sagittal l spine</td> <td>3mm x 3mm</td> <td>I70h very sharp</td> <td>pac</td> </tr> <tr> <td>5</td> <td>thin l spine</td> <td>.75mm x .7mm</td> <td>I70h very sharp</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	axial l spine tissue	3mm x 3mm	I40f medium	pac	2	axial l spine bone	3mm x 3mm	I70h very sharp	pac	3	coronal l spine	3mm x 3mm	I70h very sharp	pac	4	sagittal l spine	3mm x 3mm	I70h very sharp	pac	5	thin l spine	.75mm x .7mm	I70h very sharp	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination																											
1	axial l spine tissue	3mm x 3mm	I40f medium	pac																											
2	axial l spine bone	3mm x 3mm	I70h very sharp	pac																											
3	coronal l spine	3mm x 3mm	I70h very sharp	pac																											
4	sagittal l spine	3mm x 3mm	I70h very sharp	pac																											
5	thin l spine	.75mm x .7mm	I70h very sharp	terarecon																											
Scan Start / End Locations DFOV	external auditory meatus mid body of T1 18cm decrease appropriately																														
IV Contrast Volume / Type / Rate	70mL Iohexol (Omnipaque 300) , 2mL/sec if prescribed																														
Scan Delay	65 seconds																														
2D/3D Technique Used	Workstream 4D mpr of 3mm x 3mm coronal and sagittal l spine series, auto-transferred to PACS.																														
Comments: Recon 5 is a thin helical volume of the l spine that is archived to the TeraRecon server.																															
Images required in PACS	Topograms, 3mm x 3mm axial l spine soft tissue, 3mm x 3mm axial l spine bone, 3mm x 3mm coronal l spine, 3mm x 3mm sagittal l spine, Patient Protocol																														