




PEDIATRIC AIRWAY MANAGEMENT FOR COVID-19/PUI





LTV 1200 Circuit
(PALL Filter
Inline at the
wye)

DO NOT USE PALL
FILTER (INCREASED
DEAD SPACE) FOR
ANY TIDAL VOLUME
<350 CC
(USE HMEF FILTER)



Trilogy Circuit
(PALL Filter Inline
prior to
exhalation valve)



Inter-surgical filter
application with the
Bag

Limited number of pediatric HMEF
available. More on order. Dead space on
peds size is 37ml, infant size is 15 ml.



Stored on H2
and in HED
trauma room
airway cart

