

## RIH – IV CONTRAST NECK CHEST SIEMENS DEFINITION AS+ PROTOCOL

Indications - mass, lymphoma, adenopathy, mets.

<b>Position/Landmark</b>	Head first or feet first-Supine 1cm superior to skull vertex																									
<b>Topogram Direction</b>	Craniocaudal / Craniocaudal																									
<b>Respiratory Phase</b>	Inspiration																									
<b>Scan Type</b>	Helical																									
<b>Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization</b>	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%; text-align: center;">Neck</td> <td style="width: 50%; text-align: center;">Chest</td> </tr> <tr> <td>Care kV 120/Care Dose4D 100/0.5sec</td> <td>Care kV 120/Care Dose4D 150/0.5sec</td> </tr> <tr> <td style="text-align: center;">.8:1 , 32.00mm</td> <td style="text-align: center;">.6:1 , 24.00mm</td> </tr> <tr> <td style="text-align: center;">3 / 7</td> <td style="text-align: center;">3 / 6</td> </tr> </table>	Neck	Chest	Care kV 120/Care Dose4D 100/0.5sec	Care kV 120/Care Dose4D 150/0.5sec	.8:1 , 32.00mm	.6:1 , 24.00mm	3 / 7	3 / 6																	
Neck	Chest																									
Care kV 120/Care Dose4D 100/0.5sec	Care kV 120/Care Dose4D 150/0.5sec																									
.8:1 , 32.00mm	.6:1 , 24.00mm																									
3 / 7	3 / 6																									
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 64 = 40mm (128 x .6mm)																									
<b>Average Tube Output</b>	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%;">First Helical: ctdi – 9.7mGy dlp – 305.6 mGy.cm</td> <td style="width: 50%;">Second Helical: ctdi – 10.0mGy dlp – 366 mGy.cm</td> </tr> </table>	First Helical: ctdi – 9.7mGy dlp – 305.6 mGy.cm	Second Helical: ctdi – 10.0mGy dlp – 366 mGy.cm																							
First Helical: ctdi – 9.7mGy dlp – 305.6 mGy.cm	Second Helical: ctdi – 10.0mGy dlp – 366 mGy.cm																									
<b>First Helical Set</b> Slice Thickness/ Spacing Algorithm Recon Destination	<table style="width: 100%; border: none;"> <thead> <tr> <th style="text-align: left;">recon</th> <th style="text-align: center;">body part</th> <th style="text-align: center;">thickness/ spacing</th> <th style="text-align: center;">algorithm</th> <th style="text-align: center;">recon destination</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;"><b>axial iv neck</b></td> <td style="text-align: center;">3mm x 3mm</td> <td style="text-align: center;">I40f medium</td> <td style="text-align: center;">pacs</td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;"><b>coronal iv neck</b></td> <td style="text-align: center;">3mm x 3mm</td> <td style="text-align: center;">I40f medium</td> <td style="text-align: center;">pacs</td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">thin neck</td> <td style="text-align: center;">.75mm x .7mm</td> <td style="text-align: center;">I40f medium</td> <td style="text-align: center;">terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	<b>axial iv neck</b>	3mm x 3mm	I40f medium	pacs	2	<b>coronal iv neck</b>	3mm x 3mm	I40f medium	pacs	3	thin neck	.75mm x .7mm	I40f medium	terarecon					
recon	body part	thickness/ spacing	algorithm	recon destination																						
1	<b>axial iv neck</b>	3mm x 3mm	I40f medium	pacs																						
2	<b>coronal iv neck</b>	3mm x 3mm	I40f medium	pacs																						
3	thin neck	.75mm x .7mm	I40f medium	terarecon																						
<b>Second Helical Set</b> Slice Thickness/ Spacing Algorithm Recon Destination	<table style="width: 100%; border: none;"> <thead> <tr> <th style="text-align: left;">recon</th> <th style="text-align: center;">body part</th> <th style="text-align: center;">thickness/ spacing</th> <th style="text-align: center;">algorithm</th> <th style="text-align: center;">recon destination</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;"><b>chest</b></td> <td style="text-align: center;">5mm x 5mm</td> <td style="text-align: center;">I40f medium</td> <td style="text-align: center;">pacs</td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;"><b>lungs</b></td> <td style="text-align: center;">5mm x 5mm</td> <td style="text-align: center;">I70f very sharp</td> <td style="text-align: center;">pacs</td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;"><b>coronal chest</b></td> <td style="text-align: center;">5mm x 5mm</td> <td style="text-align: center;">I40f medium</td> <td style="text-align: center;">pacs</td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">thin chest</td> <td style="text-align: center;">.75mm x .7mm</td> <td style="text-align: center;">I40f medium</td> <td style="text-align: center;">terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	<b>chest</b>	5mm x 5mm	I40f medium	pacs	2	<b>lungs</b>	5mm x 5mm	I70f very sharp	pacs	3	<b>coronal chest</b>	5mm x 5mm	I40f medium	pacs	4	thin chest	.75mm x .7mm	I40f medium	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination																						
1	<b>chest</b>	5mm x 5mm	I40f medium	pacs																						
2	<b>lungs</b>	5mm x 5mm	I70f very sharp	pacs																						
3	<b>coronal chest</b>	5mm x 5mm	I40f medium	pacs																						
4	thin chest	.75mm x .7mm	I40f medium	terarecon																						
<b>Scan Start / End Locations</b>  <b>DFOV</b>	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%; text-align: center;">neck external auditory meatus aortic arch 18cm</td> <td style="width: 50%; text-align: center;">chest 1cm superior to lung apices 1cm inferior to adrenal glands 38cm</td> </tr> <tr> <td colspan="2" style="text-align: center;">decrease appropriately</td> </tr> </table>	neck external auditory meatus aortic arch 18cm	chest 1cm superior to lung apices 1cm inferior to adrenal glands 38cm	decrease appropriately																						
neck external auditory meatus aortic arch 18cm	chest 1cm superior to lung apices 1cm inferior to adrenal glands 38cm																									
decrease appropriately																										
<b>IV Contrast Volume / Type / Rate</b>	90mL Iohexol (Omnipaque 300) , 3mL/sec 55 second scan delay																									
<b>Scan Delay</b>	30 seconds																									
<b>2D/3D Technique Used</b>	<p>Workstream 4D mpr of 3mm x 3mm <b>coronal neck</b> series, auto-transferred to PACS.</p> <p>Workstream 4D mpr of 5mm x 5mm <b>coronal chest</b> series, auto-transferred to PACS.</p>																									
<b>Comments:</b> Recon 3 is a thin helical volume of the neck that is archived to the TeraRecon server. Recon 4 is a thin helical volume of the chest that is archived to the TeraRecon server.																										
<b>Images required in PACS</b>	Topograms, 3mm x 3mm neck, 3mm x 3mm coronal neck, 5mm x 5mm axial chest, 5mm x 5mm coronal chest, 5mm x 5mm axial lungs, Patient Protocol																									